




Only Touch What Matters

HR HealthCare offers a comprehensive range of urinary catheters, each uniquely designed to support specific clinical and patient needs. **EZ-Protect™ Closed System Intermittent Catheter** is the classic choice for no-touch, self-contained catheterization with the added versatility of two tip options.

Designed for Protection







The EZ-Protect™ Closed System Intermittent Catheter features a 1500 ml collection bag and pre-lubricated introducer tip to support hygienic catheterization. The product line is wide-ranging and includes both straight and coude tips in non-kits, standard-size kits, and pocket-size kits configurations as well as soft options.

- The soft silicone, **pre-lubricated introducer tip** helps reduce the risk of urinary tract infections by bypassing bacteria in the distal urethra.
- The catheter is self-contained within the collection bag – providing touchless use to protect the sterility of the catheter.



- **Heat-polished eyelets** for smooth insertion and withdrawal.
- Not made with natural rubber latex, DEHP, BPA, or other phthalates.

Ordering Information

EZ-Protect					Quantity: Case of 100 16" Unisex	Kit Components Conveniently packaged in order of use	
Tip	FR Size	Standard-Size Non-Kit	Standard-Size Kit [†]	Pocket-Size Kit			
straight tip (A4353)*	○ 12 Fr	20112	22112 BZK Swab	12112 BZK Swab	     	1500 ml Collection Bag Kits: EZ-Protect Closed System Catheter, Vinyl Gloves, Underpad, PVP Swabsticks or BZK Swab (as noted in chart), Gauze, and Privacy/Disposal Bag [†] Straight Tip Standard Size Kits do not include Privacy/Disposal Bag *HCPCS Codes Disclaimer: The reimbursement information provided by HR HealthCare. HR HealthCare is for convenience and general guidance only. HR HealthCare does not guarantee coverage or payment for any products.	
	● 14 Fr	20114	21114 (PVP Swabsticks) 22114 (BZK Swab)	12114 BZK Swab			
	● 16 Fr	20116	-	12116 BZK Swab			
coude tip (A4353)*	○ 12 Fr	20612	22612 BZK Swab	-			
	○ 12 Fr - Soft	-	22712 BZK Swab	-			
	● 14 Fr	20614	22614 BZK Swab	-			
	● 14 Fr - Soft	-	22714 BZK Swab	-			
	● 16 Fr	-	22616 BZK Swab	-			
	● 16 Fr - Soft	-	22716 BZK Swab	-			



To Learn More or Request Samples

Contact HR HealthCare:

hrhealthcare.com | 800.985.9048 | customercare@hrhealthcare.com

HR HealthCare, the HR HealthCare logo, and EZ-Protect are trademarks of HR HealthCare. ©2023-2024 HR HealthCare. LBL-0080 Rev. 1 (04/2024)