

At its most fundamental level, **self-determination** is the right of people to shape their own lives and live well.




## Self-Determination. Don't Hold Back.

Stretching the limits of what may have been previously impossible, the **EZ-Gripper™ Closed System Intermittent Catheter** and **Eagle Board™ Male Catheter Assistance Device** were engineered to support limited dexterity user needs and facilitate one-handed catheterization. The result is a unique solution to empower confidence, foster self-determination, and forge independence.


# Discover the Power Within

The EZ-Gripper and Eagle Board provide innovative features to champion and assist with independent self-catheterization. Discover the unparalleled combination of features and challenge the possibility of establishing a self-supporting bladder management plan for greater flexibility and freedom.



Use the EZ-Gripper alone or with the Eagle Board, making one-hand catheterization possible for users with limited hand dexterity.

- Wings rest on the thighs providing stabilization and angling to support urine flow.
- Penile lever secures the penis to the board to hold it in place.
- Housing clip secures and holds the catheter bag during use.



Designed exclusively for use with EZ-Gripper

## Ordering Information

EZ-Gripper			Quantity: Case of 50 16"   18" Unisex
Length	FR Size	Standard-Size Non-Kit	Standard-Size Kit
16"	○ 12	50112	52112 BZK Swab
	● 14	50114	52114 BZK Swab
	● 16	-	52116 BZK Swab
18"	○ 12	-	52812 BZK Swab

### Kit Components

Conveniently packaged in order of use



1500 ml Collection Bag  
Kits: EZ-Gripper Closed System Catheter, Vinyl Gloves, Underpad, BZK Swab, Gauze, and Privacy/Disposal Bag.

\*A4353 - HCPCS Codes Disclaimer: The reimbursement information provided by HR HealthCare is for convenience and general guidance only. HR HealthCare does not guarantee coverage or payment for any products.

### Eagle Board

REF# 100MTG

Includes: Eagle Board with Penile Sizing Guide, Foam Pads, Carry Bag, Allen Wrench, and Instructions for Use.

Not made with natural rubber latex, DEHP, BPA, or other phthalates.



To Learn More or Request Samples

Contact HR HealthCare:

hrhealthcare.com | 800.985.9048 | customercare@hrhealthcare.com

HR HealthCare, the HR HealthCare logo, EZ-Gripper, EZ Advancer, and Eagle Board are trademarks or registered trademarks of HR HealthCare. ©2023-2024 HR HealthCare. All rights reserved. LBL-0082 Rev. 2 (04/2024)

