



THE CLASSIC NO-TOUCH CLOSED SYSTEM CATHETER

For those who prefer the classic, closed system design, MTG offers our high quality MTG Jiffy Cath[®] no-touch closed system catheter. Designed to help reduce the risk of recurrent urinary tract infections with the all-important introducer tip, the MTG Jiffy Cath[®] has soft, fire-polished eyelets for a smooth, comfortable catheterization experience.

Available in 12-16 FR sizes and 16" length.

MTG Jiffy Cath[®]

No-Touch Closed System Catheter

FEATURES & BENEFITS

Designed to Reduce Risk of Recurrent UTIs

The no-touch closed system catheter and pre-lubricated introducer tip are designed to reduce the risk of recurrent urinary tract infections.

Collection Bag Protects Sterility

The self-contained 1500ml collection bag protects the sterility of the catheter tube and collects urine while measuring output.

Available in Straight Tip and Coudé Tip

For those who have conditions that can make inserting a straight tip catheter difficult, we offer MTG Jiffy Cath[®] coudé tip closed system catheters. Available in soft or firm catheter tube options. *Soft coudé is available in standard kits only.*

Self-Cath Anywhere with Mini-Pak Size Kits

Our Mini-Pak size kits are the perfect traveling option and fit discreetly in your pocket or purse. Mini-Pak size kits include a privacy bag for easy and discreet disposal.

Made with Medical Grade Quality Materials

All MTG catheters use medical grade quality materials and are not made with natural rubber latex, BPA or DEHP.





MTG catheter kits include everything you need to self-cath. Full kits contain one sterile MTG Jiffy Cath® closed system catheter, pair of latex-free gloves, underpad, gauze pad, Povidone-iodine (PVP) or Benzalkonium chloride (BZK) clear, odor-free skin cleansing prep wipe and privacy bag.

MTG Jiffy Cath® Straight Tip No-Touch Closed System Catheters

Single Closed System Catheters

Item #	FR Size	Description	Antiseptic Wipe
20112	12	MTG Jiffy Cath® No-Touch Closed System Catheter, 12 Fr, 16"	NA
20114	14	MTG Jiffy Cath® No-Touch Closed System Catheter, 14 Fr, 16"	NA
20116	16	MTG Jiffy Cath® No-Touch Closed System Catheter, 16 Fr, 16"	NA

Closed System Catheter Kits - Standard Size

Item #	FR Size	Description	Antiseptic Wipe
22112*	12	MTG Jiffy Cath® No-Touch Closed System Catheter Kit, 12 Fr, 16"	BZK Wipe
21114*	14	MTG Jiffy Cath® No-Touch Closed System Catheter Kit, 14 Fr, 16"	PVP Swabs
22114*	14	MTG Jiffy Cath® No-Touch Closed System Catheter Kit, 14 Fr, 16"	BZK Wipe
22116*	16	MTG Jiffy Cath® No-Touch Closed System Catheter Kit, 16 Fr, 16"	BZK Wipe

Closed System Catheter Kits - Mini-Pak Size

Item #	FR Size	Description	Antiseptic Wipe
12112	12	MTG Jiffy Cath® No-Touch Closed System Catheter Kit, Mini-Pak Size, 12 Fr, 16"	BZK Wipe
12114	14	MTG Jiffy Cath® No-Touch Closed System Catheter Kit, Mini-Pak Size, 14 Fr, 16"	BZK Wipe
12116	16	MTG Jiffy Cath® No-Touch Closed System Catheter Kit, Mini-Pak Size, 16 Fr, 16"	BZK Wipe

MTG Jiffy Cath® Coudé Tip No Touch Closed System Catheters

Single Closed System Catheters

Item #	FR Size	Description	Antiseptic Wipe
20612	12	MTG Jiffy Cath® No-Touch Closed System Catheter, Coudé, 12 Fr, 16"	NA
20614	14	MTG Jiffy Cath® No-Touch Closed System Catheter, Coudé, 14 Fr, 16"	NA

Closed System Catheter Kits - Standard Size

Item #	FR Size	Description	Antiseptic Wipe
22612	12	MTG Jiffy Cath® No-Touch Closed System Catheter Kit, Coudé, 12 Fr, 16"	BZK Wipe
22614	14	MTG Jiffy Cath® No-Touch Closed System Catheter Kit, Coudé, 14 Fr, 16"	BZK Wipe
22616	16	MTG Jiffy Cath® No-Touch Closed System Catheter Kit, Coudé, 16 Fr, 16"	BZK Wipe

Soft Coudé Closed System Catheter Kits

Item #	FR Size	Description	Antiseptic Wipe
22712	12	MTG Jiffy Cath® No-Touch Closed System Catheter Kit, Soft Coudé, 12 Fr, 16"	BZK Wipe
22714	14	MTG Jiffy Cath® No-Touch Closed System Catheter Kit, Soft Coudé, 14 Fr, 16"	BZK Wipe
22716	16	MTG Jiffy Cath® No-Touch Closed System Catheter Kit, Soft Coudé, 16 Fr, 16"	BZK Wipe

*Indicates items that do not include privacy bags.

HCPCS: A4353, Packaging: 100/Case